

Liver Assist


Unlock the power of flexibility



Hypothermic Oxygenated Perfusion raising HOPE for donated livers

Results from clinical trials demonstrate that hypothermic machine perfusion of livers leads to:

- Improved 5-year patient and graft survival in HOPE treated DCD livers vs. untreated DCD livers¹
- Lower biliary complications in dHOPE treated DCD livers compared to untreated DCD livers^{1,2,3}
- Less ischemia reperfusion injury after oxygenated perfusion in DCD livers²
- Significant shorter ICU and hospital stay in ECD grafts after HOPE⁴
- End-ischemic HOPE reduces early allograft injury and improves transplant outcomes in ECD-DBD liver transplantation⁵



Normothermic ex vivo liver perfusion for increased utilization

Normothermic liver perfusion leads to increased availability by viability assessment of liver grafts. Published data suggests that normothermic ex vivo liver perfusion leads to:

- 20% increase in liver transplantation combining hypothermic and normothermic machine perfusion⁶
- Safe use of initially rejected donated livers^{7,8}

Liver Assist

Unlock the power of flexibility

With more than 10 years of clinical application, XVIVO's Liver Assist is the most used device for ex vivo perfusion of livers. The two separately controlled pump units provide oxygenated perfusion with near physiologic settings with pulsatile perfusion of the hepatic artery and continuous flow through the portal vein.

XVIVO's Liver Assist automatically adjusts the flow to the natural resistance of the organ. The heater/cooler unit enables perfusion at every temperature between hypothermic and normothermic.

With adjustable settings, XVIVO's Liver Assist allows for the clinician's choice of protocol, including HOPE, dHOPE, COR and NMP.



Oxygenated



Pressure controlled
pulsatile (60 bpm) and
continuous perfusion



Temperature
controlled perfusion
12°C - 37°C



Functional viability
assessment during NMP
(bile and perfusate analysis)



Up to 24H cold
perfusion to
support logistics

Oxygenation

Oxygenation and/or gas exchange via two separate hollow fiber oxygenators.

Pulsatility

The hepatic artery is perfused with a 60 bpm pulse generated by a centrifugal pump.

Flexible temperature

The only device that allows perfusion within a flexible temperature range of 12-37 °C.

Ease of use

Simple user interface for easy operation. Provides an ergonomic working height.

Filter

Perfusate is filtered before it enters the liver.

Organ chamber


Easy to access the organ during perfusion, while ensuring a safe, sterile and humid environment. The reservoir has an integrated bile drain.

Monitoring & sampling

Integrated sensors monitor real time flows, pressures and temperatures. Incorporated sampling port allows for analysis of perfusate and bile (such as pH, lactate, glucose, enzymes and others).

Continuous

Mimicking physiological circumstances the portal vein is perfused with a continuous flow generated by a centrifugal pump.



The DHOPE-DCD trial: new hope for DCD liver perfusion³


The study aimed to compare hypothermic oxygenated machine perfusion (2h, end-ischemic) to static cold preservation of livers donated after circulatory death (DCD). The study was: randomized, controlled, multicenter, n=156 (78+78).

The study showed that:

- Nonanastomotic strictures occurred in 6% of the patients in the dHOPE group and in 18% in the control group (risk ratio 0.36, P=0.03).
- The cumulative number of treatments for nonanastomotic biliary strictures was lower by a factor of nearly 4 in the dHOPE group compared to the control group.

According to the authors, the prevention of post-transplant cholangiopathy may increase the acceptance for DCD livers and make the use of machine perfusion cost-effective.

Biliary complications significantly reduced with hypothermic oxygenated machine perfusion of DCD livers³



Cumulative Incidence of Symptomatic Nonanastomotic Biliary Strictures: adapted from³ Hazard ratio, 0,32 (95% CI, 0.11 - 0.89); P=0.03

- HOPE : Hypothermic Oxygenated Perfusion
- NRP: Normothermic Regional Perfusion
- ECD: Extended Criteria Donor
- DCD: Donation after Circulatory Death

Nobody should die waiting for a new organ

Founded in 1998, XVIVO is the only MedTech company dedicated to extending the life of all major organs – so transplant teams around the world can save more lives.

Our solutions allow leading clinicians and researchers to push the boundaries of organ transplantation.

Scan this QR-code to get in touch with us!



Not all XVIVO products are approved in all markets.

References:

1. Schlegel A, et al. Outcomes of DCD liver transplantation using organs treated by hypothermic oxygenated perfusion before implantation. *J Hepatol*. 2019 Jan;70(1):50-57.
2. Dutkowski P, et al. First Comparison of Hypothermic Oxygenated Perfusion Versus Static Cold Storage of Human Donation After Cardiac Death Liver Transplants: An International-matched Case Analysis. *Ann Surg*. 2015 Nov;262(5):764-70; discussion 770-1.
3. van Rijn R, et al. Hypothermic Machine Perfusion in Liver Transplantation – A Randomized Trial. *N Engl J Med*. 2021; 384:1391-1401.
4. Rayar M, et al. Hypothermic Oxygenated Perfusion Improves Extended Criteria Donor Liver Graft Function and Reduces Duration of Hospitalization Without Extra Cost: The PERPHO Study. *Liver Transpl*. 2021 Feb;27(3):349-362.
5. Czigan Z, et al. Hypothermic Oxygenated Machine Perfusion (HOPE) Reduces Early Allograft Injury and Improves Post-Transplant Outcomes in Extended Criteria Donation (ECD) Liver Transplantation from Donation After Brain Death (DBD): Results from a Multicenter Randomized Controlled Trial (HOPE ECD-DBD). *Ann Surg*. 2021 Jul 29.
6. van Leeuwen OB, et al. Transplantation of High-risk Donor Livers After Ex Situ Resuscitation and Assessment Using Combined Hypo- and Normothermic Machine Perfusion: A Prospective Clinical Trial. *Ann Surg*. 2019 Nov;270(5):906-914.
7. Watson C, et al. Normothermic Perfusion in the Assessment and Preservation of Declined Livers Before Transplantation: Hyperoxia and Vasoplegia—Important Lessons From the First 12 Cases. *Transplantation*. 2017 May; 101(5): 1084-1098.
8. Mergental H, et al. Transplantation of Declined Liver Allografts Following Normothermic Ex-Situ Evaluation. *Am J Transplant*. 2016 Nov;16(11):3235-3245.
9. Watson C, et al. Observations on the ex situ perfusion of livers for transplantation. *Am J Transplant*. 2018 Aug;18(8):2005-2020.

XVIVO BV, Bornholmstraat 84, 9723 AZ Groningen, Netherlands
Tel +31 50 313 1905 | xvivogroup.com

XVIVO trademarks, names and logos are registered trademarks of XVIVO Perfusion AB.
© 2023 XVIVO Perfusion AB. All rights reserved.

TD-11 Brochure LiA.09

XVIVO

XVIVO